



**RESTORATIVE ARTS**  
DENTAL LABORATORY

33B Gamecock Ave  
Charleston, SC 29407  
800-333-6527 843-571-6656  
rad@restorativearts.com

The compatible components listed fit directly to the implant. For example, the UCLA components for the Branemark RP may also be used with the 3i miniimplant & standard 3.75 & 4.0, Threadloc, etc. The external configuration of the compatible implants may differ, but the fitting surface and internal threads are compatible.

Implants	Restoratively Compatible To	Implants	Restoratively Compatible To
<b>Branemark</b>	<b>Branemark NP</b> 3.3 external hex - 3.5 Platform <b>Branemark RP</b> 3.7 & 4.0 external hex - 4.1 Platform	<b>Restorative SD</b> threaded or cylindrical ex hex 3.3 or 3.4 <b>Sustain SD</b> threaded or cylindrical ex hex 3.3 or 3.4 <b>Restore RD</b> threaded or cylindrical ex hex 3.75 or 4.0 <b>Sustain RD</b> threaded or cylindrical ex hex 3.75, 4.0, 4.2 <b>Stage-1</b> internal locking Morse taper <b>Sustain</b> internal bevel cylindrical 3 diameters	<b>Restoratively Compatible To</b> Sustain SD 3.4 Restore SD 3.3 Branemark, Regular Platform 4.1 Branemark, Regular Platform 4.1 ITI none
<b>Branemark</b>	<b>Branemark MK III &amp; IV</b> 3.3 Implanted NP <b>Branemark MK III &amp; IV</b> 3.75 & 4.0 NobelSpeedy Groovy NobelSpeedy Groovy Shorty 3i Miniimplant & Standard 3.75 & 4.0 Implanted RP <b>IMZ Hex Head</b> 3.3 & 4.0 Innova Entegra 3.25 & 4.0 Osseo Implant 3.75 Restore RD & Sustain RD Steri-Oss Hex-Loc 3.8 & 4.5 Screw-Vent 3.75 Taper-Lock Threadloc	<b>Core-Vent SD</b> 3.5 <b>Cone-Vent MD</b> 4.5 <b>Screw-Vent WD</b> 5.7 - 6.0 Platform <b>Taper-Lock</b> external hex - 4.1 Platform <b>Swiss Plus</b>	<b>Restoratively Compatible To</b> IMZ 3.5, SV 3.7 & 3.75, BV 3.5, ADI 3.5 MV 4.25 & 4.7, SV 4.7, BV 4.5, ADI 4.5 none Branemark, Regular Platform 4.1 ITI Syn-Octa 4.8, Blue Sky Bio
<b>Branemark</b>	<b>WP</b> 5.0 & 5.5 external hex - 5.1 Platform	<b>LifeCore</b> <b>3i TG Osseofit</b>	<b>Restoratively Compatible To</b> Branemark MK III & IV 5.0 Implanted WP
<b>IMZ</b>	<b>IMZ Original</b> 3.3 & 4.0 non-hex <b>IMZ Hex Head</b> 3.3 & 4.0 external hex	<b>Steri-Oss Original</b> 3.5, 3.8 and 4.0 non-hex flat-top <b>Steri-Oss Hex-Loc</b> 3.25 <b>Steri-Oss Hex-Loc</b> 3.8 & 4.5 <b>Steri-Oss Replace</b> 3.5 & 4.3 external hex <b>Steri-Oss Replace</b> 5.0 external hex <b>Steri-Oss Select</b> internal connection	<b>Restoratively Compatible To</b> none O-Burment / Anchor 4.0 Integral 4.0 none Branemark, Regular Platform 4.1 none Steri-Oss HexLoc 5.0 NobelReplace Tapered Groovy NobelSpeedy Replace NobelReplace Straight
<b>ITI</b>	<b>ITI Small Diameter</b> external hex <b>ITI Standard</b> 4.8 internal Morse taper	<b>Paragon SwissPlus</b> , Blue Sky Bio	<b>Restoratively Compatible To</b> none LifeCore Stage-1 3i TG Osseofit
<b>Integral</b>	<b>ITI Syn-Octa</b> 4.8 internal octagonal Morse taper <b>Integral SD</b> 3.25 non-hex flat-top <b>Integral</b> 4.0 non-hex flat-top <b>Omniloc</b> 3.25 & 4.0 <b>Threadloc</b> 3.75 external hex <b>Calcitek Spline SD, RD &amp; WD</b>	<b>O-Burment</b> 3.25 O-Burment / Anchor 4.0 Steri-Oss Original 4.0 non-hex none Branemark, Regular Platform 4.1 none	<b>Restoratively Compatible To</b> none O-Burment / Anchor 4.0 Integral 4.0 none Branemark, Regular Platform 4.1 none Steri-Oss HexLoc 5.0 NobelReplace Tapered Groovy NobelSpeedy Replace NobelReplace Straight
<b>3i-Implant Innovations</b>	<b>3i Microimplant</b> 3.25 external hex 3.4 platform <b>3i Miniimplant</b> 3.25 external hex 4.1 platform <b>3i Standard</b> 3.75 external hex 4.1 platform <b>3i Wide Diameter</b> 5.0 external hex 5.0 platform <b>3i Wide Diameter</b> 6.0 external hex 6.0 platform <b>3i Osseofit TG</b> internal Morse taper	<b>ITI Standard</b> 4.8 internal taper	<b>Restoratively Compatible To</b> none Branemark, Regular Platform 4.1 Branemark, Regular Platform 4.1 none none ITI Standard 4.8 internal taper

## Compatible Implant Restorative Components